

MOBILETT XP Digital

SP

Maintenance Protocol  
System

incl. DHHS

Customer:

Address:

Department:

Room:

Contact person:

Telephone:

Cust. specific no.:

Cust. no.:

Date.:

The instructions SPR8-230.831.30.03.02 are required for  
this protocol

© Siemens AG 2004  
The reproduction, transmission or use  
of this document or its contents is not  
permitted without express written  
authority. Offenders will be liable for  
damages. All rights, including rights  
created by patent grant or registration  
of a utility model or design, are  
reserved.

<b>SIEMENS Office:</b>	
<b>Address:</b>	
<b>Region:</b>	
<b>Country:</b>	
<b>Contact person:</b>	
<b>Tel.:</b>	
<b>CSE in charge:</b>	
<b>Tel.:</b>	

**Remarks Regarding the Protocol:**

The protocol is valid as proof of quality for **one** check that must be performed on the system / component in one year.

The check must be performed in the specified intervals.

The results of the check are entered in this protocol.

The chapter numbers in front of the checkpoints indicate the corresponding chapters in the particular instructions (see cover page).

The protocol must be completely filled out by the Customer Service Engineer, i.e.:

- All boxes must be filled out. If a box does not apply to the system or if no entry needs to be made, check the "n.a." box.
- Enter the customer number (Cust. No. :) and the date of the check in the header of each page so that each page can be allocated to a customer and to a check date.
- If there are complaints, the IVKs for the component about which a complaint has been made as well as the type of complaint must be entered in the "Open Points" table provided for this. Correction of these open points also must be documented in this table with the date and a signature. If there are no open points, check "No" and document this with the date and a signature.
- If movable components (also test phantoms that are part of the system) that can be used in different systems are used for the check, they must be entered in the "Movable Components" table provided for this.
- The measurement values for the measurements that must be performed during the check must also be entered in the open spaces / tables provided for them.
- After completing the check, Page 3 of this protocol must be filled out and signed.

**Further Processing and Archiving of the Protocol**

The protocol is a document and thus must be archived. After completing the test, it must be filed in the corresponding register in the "System Owner Manual" binder. If needed, a copy can be handed to the customer.

<b>System:</b>	
<b>Serial No.:</b>	
<b>Software Version:</b>	
<b>Number of the Service Contract:</b>	
<b>Type of Maintenance:</b>	

**Evaluating the Condition of the System / Component**

The system has no deficiencies. The image quality test resulted in no differences from required reference values.	
The system / component has slight deficiencies that have no affect on continued operation of the system. However they should be corrected preventively. The image quality test resulted in no differences from required reference values.	
The system / component has serious deficiencies. For safety reasons, continued operation of the system is permitted only after successfully correcting the deficiencies.	

**After completing all work steps, an evaluation was performed.**

Signature: \_\_\_\_\_

Date:

Name:

The operator or a person assigned for this has taken note of this evaluation.  
(if national regulations require this)

Signature: \_\_\_\_\_

Date:

Name:

## Explanation of Abbreviations in the Protocol

Abbrev.	Explanation	Abbrev.	Explanation
SI	Safety Inspection	PMF	Preventive Maintenance, Operating Value Check, Function Check
SIE	Electrical Safety Inspection	Q	System Quality, Image Quality
SIM	Mechanical Safety Inspection	QIQ	Image Quality
PM	Preventive Maintenance	QSQ	System Quality Check
PMP	Periodic Preventive Maintenance	SW	Software Maintenance
PMA	Preventive Maintenance Adjustments	CSE	Customer Service Engineer

## Additional activities performed

Only activities that are not described in the instructions for the system / component need to be listed.

	Date:		
Additional activities performed:	OK	not OK	n.a.

## Open Points:

Yes:                      No:                      Signature: \_\_\_\_\_

                                    Date:                      Name: \_\_\_\_\_

If "Yes", enter the component with the IVK and the open point (only the number) in the table. After completing maintenance, record this in the table.

IVK	Component	Open Points	Completed	
			Date	Signature

**Measuring Devices queried electronically:**

Yes:                      No:                      Signature: \_\_\_\_\_  
Date:                      Name:

If the measurement devices are queried electronically, for example with a Scout Mobile Device, entry of the measuring devices in the table can be skipped.

Measuring Devices	Type	Serial No.	Date Used	Next Calibration Due

**Movable Components:**

Yes:                      No:                      Signature: \_\_\_\_\_  
Date:                      Name:

If "Yes", enter the movable component with which the check was performed along with the Serial No. in the table.

Movable components (also test phantoms that are part of the system) are parts that can be used on different systems).

Component	Serial No.

- 1 General**
- 2 General maintenance information**
- 3 Inspection and maintenance**
  - 3.1 Visual check**
    - SIM Signs
    - SI Customer documentation
  - 3.1.1 Damage**
    - SIM Covers
    - SIM Detector and detector holder
    - SIM Control panel and display panel
    - SIM Release cable
    - SIM Single tank
    - SIM Collimator
    - SIM DAP (optional)
    - SIM Remote control (optional)
  - 3.2 Mechanical inspection**
    - 3.2.1 Back wheels and support rollers**
      - 3.2.1.1 Back wheels**
        - SIM Brake pad
        - SIM Drive belt
        - SIM Secure attachment
        - SIM Smooth rotation
      - 3.2.1.2 Support rollers**
        - SIM Secure attachment
        - SIM Smooth rotation
    - 3.2.2 Front transport wheels**
      - SIM Secure attachment
      - SIM Smooth rotation
    - 3.2.3 Brakes**
      - SIM Uniformity
      - SIM Braking force
    - 3.2.4 Support arm transport lock**
      - SIM Checking transport lock
      - SIM Locking mechanism
    - 3.2.5 Handles**
      - SIM Attachment

**3.2.6 Collimator adjustment knobs**

SIM Attachment

**3.2.7 Arm system and single tank****3.2.7.1 Moving the arm system**

SIM Moving the arm system

SIM Moving the single tank holder

SIM Moving the single tank

**3.2.7.2 Arm system attachment**

SIM Screws and lock nuts

**3.2.7.3 Arm connector**

SIM Screws and nuts

**3.2.7.4 Adjusting screw and cantilever**

SIM Lock nuts

SIM Attaching the base

**3.2.7.5 Cable harness for the arm system**

SIM Undamaged

**3.2.8 Single tank holder**

SIM Holder screw connections

SIM Single-tank screw connections

**3.2.9 Power cable**

PMP Cable winch replacement

Start-up date/date of last replacement:

. .  
DD MM YYYY

SIM Checking the power cable

SIM General cable winch information

PMP Cleaning

**3.2.9.1 Power plug**

SIM Checking the power plug

**3.2.10 Lubrication**

PMP Pull-bar slide bushing

**3.3 Function inspection****3.3.1 Operating data****3.3.1.1 Reading out the operating information**

PMF Exposure counter

Measured value:

PMF Error memory

**3.3.2 Displaying the control panel**

PMF Mains display

PMF kV/mAs default values

SIE kV/mAs segment displays

PMF Ready indicator

SIE kV/mAs displays according to operating instructions

**3.3.3 Checking the radiation indicator**

SIE Radiation indicator

SIE Ready indicator

SIE Acoustic signal

**3.3.4 Manual termination of exposure**

SIE "ERR 39" display

SIE Acoustic signals

**3.4 Collimator****3.4.1 Lamp replacement**

PMP Annual replacement of the collimator lamp

Start-up date/date of last replacement: . .

DD MM YYYY

PMF Light localizer function

**3.4.2 Checking the illuminance**

PMF Illuminance

Measured value:

**3.4.3 Light field/radiation field**

QSQ Deviation ((A + C) / SID)

Measured value:

QSQ Deviation ((B + D) / SID)

Measured value:

**3.5 Battery and motor drive inspection****3.5.1 Batteries**

PMF Charging performance

Start-up date/date of last replacement: . .

DD MM YYYY



**3.5.2 Motor drive**

PMF Forward/backward

PMF Slow/fast

**3.6 Options****3.6.1 DAP measuring system**

SIE Function of the DAP measuring system

**3.6.2 Remote control**

PMP Battery replacement for remote control (yearly)

SIE Remote control function

**3.7 Checking the kV/mAs exposure parameters****3.7.1 kV accuracy**

PMF kV accuracy 52 kV, 50 mAs

Measured value:

PMF kV accuracy 81 kV, 20 mAs

Measured value:

PMF kV accuracy 133 kV, 12.5 mAs

Measured value:

**3.7.2 mAs accuracy**

PMF kV accuracy 40 kV, 5 mAs

Measured value:

PMF kV accuracy 81 kV, 2 mAs

Measured value:

PMF kV accuracy 133 kV, 10 mAs

Measured value:

**3.8 Reproducibility test (USA only)**

QSQ Determine coefficient of variation C.

Measured value:

**3.9 Checking the image quality****3.9.1 Dose measurement**

QIQ Dose measurement

**3.9.2 Resolution**

QIQ Resolution

**3.9.3 Contrast**

QIQ Contrast

**3.9.4 Hardcopy**

QIQ Hardcopy

**3.10 Protective conductor test**

SIE Protective conductor test

**3.11 Leakage current measurement**

SIE Leakage current measurement

**3.12 Patient leakage current measurement**

SIE Patient leakage current measurement

**3.13 Cleaning**

PMP Cleaning